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


## Mepitel® Ag

Antimicrobial soft silicone wound contact layer

**Read only**
**Caution: Federal [US] Law restricts this device to sale by or on the order of a physician [or properly licensed practitioner].**

### Product description

Mepitel Ag is a soft silicone wound contact layer that allows exudate to pass vertically into a secondary dressing.

Mepitel Ag makes it possible to change only the secondary absorbent dressing. It maintains its structural integrity and can be left in place for up to eight (8) days depending on the wound condition and surrounding skin (or as indicated by accepted clinical practice).

In vitro testing demonstrates that Mepitel Ag gives a >4 log reduction, on the dressing, of the following Gram positive bacteria, Gram negative bacteria and yeast: *Enterococcus faecalis* [VRE], *Staphylococcus aureus* (MRSA), *Staphylococcus epidermidis*, *Acinetobacter baumannii*, *Enterobacter cloacae*, *Pseudomonas aeruginosa*, *Candida albicans*, *Candida guilliermondi*, *Candida lusitiane*. The dressing sustains antimicrobial activity for up to eight (8) days in in vitro studies. By reducing the number of micro-organisms, Mepitel Ag may also reduce odour.

Mepitel Ag has shown to have a 4 log microbial reduction up to eight (8) days when tested in vitro with a secondary dressing.

Mepitel Ag maintains a moist wound environment in combination with a secondary dressing.

Mepitel Ag can be used under compression bandaging.

Mepitel Ag consists of a polyamide net with Safetac adhesive layer, containing silver sulfate and cellulose compound, on both sides.

Safetac® is a proprietary adhesive technology developed by Mölnlycke.

Dressing material content: Polydimethylsiloxane, Polyamide, Silver sulfate, Sodium carboxymethylcellulose.

#### Indications for use

Mepitel Ag is intended for the management of a wide range of exuding wounds such as skin tears, skin abrasions, sutured/surgical wounds, partial thickness burns, partial and full thickness grafts, lacerations, diabetic ulcers, venous ulcers, and arterial ulcers. Silver sulfate is added to the dressing as a preservative to inhibit or reduce microbial growth on the dressing.

#### Precautions

- Do not use Mepitel Ag on patient and/or user with known hypersensitivity to the ingoing materials/components of the product.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Frequent or prolonged use of this product may result in transient or permanent discoloration of skin.
- The silver sulfate component in Mepitel Ag is intended to act as a preservative for the dressing material. It is not intended to treat infection. In the event of clinical infection Mepitel Ag does not replace the need for systemic therapy or other adequate infection treatment.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.
- Avoid placing pressure upon the dressing when used on partial thickness burns with high risk of rapid granulation or after facial resurfacing.
- When used on bleeding wounds or wounds with high viscosity exudate, Mepitel Ag should be covered with a moist absorbent dressing.
- When Mepitel Ag is used for the fixation of skin grafts, the dressing should not be changed before the fifth day post application.

#### Instructions for use

Mepitel Ag should be used by or under the supervision of a qualified health care professional.

Note that local hygiene procedures should be followed prior to and following the dressing change.

- Cleanse the wound in accordance with normal procedures and dry the surrounding skin thoroughly.
- Select appropriate dressing size. For best clinical outcomes, ensure that the product covers the wound bed and overlaps the wound edges. Mepitel Ag can be cut to size to suit various wound shapes and locations.
- While holding the larger protective film, remove the smaller one. If needed, moisten gloves to avoid adherence to Mepitel Ag.

4. Apply Mepitel Ag over the wound and remove the remaining protective film. Smooth Mepitel Ag in place onto the surrounding skin, ensuring a good seal. If more than one piece of Mepitel Ag is required, overlap the dressing edges, ensuring that the holes are not blocked.

5. Apply a secondary dressing on top of Mepitel Ag. In contoured or jointed areas (e.g. under arm, under breast, inner elbow, groin, deep wounds), ensure sufficient padding is applied to keep the Mepitel Ag held flat against the surface of the wound.

6. Apply a suitable fixation.

Mepitel Ag

Duration of treatment is determined by the physician and depends on the wound type and healing conditions. Each Mepitel Ag may be left in place for up to eight (8) days depending on the patient, condition of the wound and surrounding skin, or as indicated by accepted clinical practice (exudates should pass freely through the dressing and the holes should not be blocked).

If saturated, the secondary absorbent dressing could be changed with the Mepitel Ag left in place.

If the condition of the wound deteriorates unexpectedly, consult a healthcare professional for the appropriate medical treatment.

The interaction of Mepitel Ag together with cleansing agents and topical treatments has not been demonstrated.

Prior to commencing radiation therapy remove Mepitel Ag if product is present in the treatment area. A new dressing can be applied following treatment.

If present during conventional X-ray imaging Mepitel Ag will appear as a grid like pattern artifact.

Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).

#### Special storage conditions and handling conditions

Mepitel Ag should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. Disposal should be handled according to local environmental procedures.

#### Other information

Mepitel Ag contains 1.8–3.0 mg/cm² silver.

If any serious incident has occurred in relation to the use of Mepitel Ag it should be reported to Mölnlycke Health Care.

Mepitel Ag

Mepitel® and Safetac® are registered trademarks of Mölnlycke Health Care AB.

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