

Mölnlycke®

Mepilex® Ag

WITH SAFETAC® TECHNOLOGY

STERILE EO



CE 0086

Manufacturer / Fabricant




 Mölnlycke Health Care AB  
Gamlestadvägen 3C, Box 13080  
SE-402 52 Göteborg, Sweden

[www.molnlycke.com](http://www.molnlycke.com)

## Mepilex® Ag

## Mepilex® Heel Ag/ Mepilex® Talon Ag

**en** For low exuding wounds

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# Mepilex® Ag with Safetac® Technology

Antimicrobial soft silicone foam dressing

## Product description

Mepilex Ag consists of:

1. a Safetac® wound contact layer
2. a flexible absorbent pad of grey polyurethane foam containing a silver compound and activated carbon
3. an outer film which is vapour permeable and waterproof

Mepilex Ag contains silver sulphate that releases silver ions to create an effective bacterial barrier and inactivates a wide range of wound related pathogens (bacteria and fungi), shown in vitro. By reducing the number of microorganisms, Mepilex Ag may also reduce odour.

Mepilex Ag has also been shown to inactivate wound related pathogens, up to 7 days in vitro.

Safetac is a unique adhesive technology that minimises pain to patients and trauma to intact skin or wounds.

## Mode of action

Mepilex Ag is a highly conformable dressing that absorbs exudate and maintains a moist wound environment.

As Mepilex Ag maintains a moist wound environment, supporting debridement, there might be an initial increase in the wound size. This is normal and to be expected.

## Intended use

Mepilex Ag is an antimicrobial soft silicone foam dressing that is designed for the management of low to moderately exuding wounds such as leg and foot ulcers, pressure ulcers and partial thickness burns.

Mepilex Ag may be used on infected wounds as part of a treatment regimen under supervision of a qualified health care professional.

Mepilex Ag can be used under compression bandaging.

## Instructions for use

Note that local hygiene procedures should be followed prior to and following the dressing change.

1. Cleanse the wound with saline solution or water according to standard clinical practice.
2. Dry the surrounding skin thoroughly.
3. Remove the release films and apply the adherent side to the wound. Do not stretch.
4. For best result, Mepilex Ag should overlap the dry surrounding skin by at least 1-2 cm for the smaller sizes (sizes up to 12.5x12.5 cm) and 5 cm for the larger sizes in order to protect the surrounding skin from maceration and excoriation and fixate the dressing securely. If required, Mepilex Ag can be cut to suit various wound shapes and locations.
5. When necessary, fixate Mepilex Ag with a bandage or other fixation.

Mepilex Ag is intended for short-term use up to 4 weeks. For long-term use, a clinical assessment by a physician is recommended.

## Frequency of change

Mepilex Ag may be left in place for several days depending on the condition of the wound and surrounding skin, or as indicated by accepted clinical practice.

A change in dressing regimen can result in an initial increased level of exudates, which temporarily may require an increased change frequency.

## Precautions

- Mepilex Ag should be used under the supervision of a qualified health care professional.
- Do not use on patients with a known sensitivity to silver.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Do not use Mepilex Ag during radiation treatment or examinations e.g. X-ray, ultrasound, diathermy or Magnetic Resonance Imaging.
- Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).

- Do not use Mepilex Ag together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- For external use only.
- Mepilex Ag may cause transient discoloration of the wound bed and surrounding skin.
- In the event of clinical infection Mepilex Ag does not replace the need for systemic therapy or other adequate infection treatment.
- The interaction of Mepilex Ag with other topical treatments has not been demonstrated.
- Other than saline solution or water, the interaction of cleansing agents in combination with Mepilex Ag has not been demonstrated.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if inner package is damaged or opened prior to use. Do not re-sterilise.
- If the product is used after the expiry date product properties cannot be ensured.

#### **Storage and disposal**

- Mepilex Ag should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. Note that any colour variances in Mepilex Ag do not affect product performance or safety.
- Disposal should be handled according to local environmental procedures.

#### **Silver content information**

- Mepilex Ag contains 1.2 mg/cm<sup>2</sup> silver

Mepilex and Safetac are registered trademarks of Mölnlycke Health Care AB.

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**CANADA** Mölnlycke Health Care Inc.,  
2010 Winston Park Drive, Suite 100,  
Oakville, Ontario, L6H 5R7  
Tel: +1 905 829 1502

**CHINA** Mölnlycke Health Care,  
Room 301, Building H,  
Beijing Phoenix Plaza, No. 5A Shuguangxili  
Chaoyang District, 100028 Beijing, P.R. China  
Tel: +86 10 5128 8571

**CZECH REPUBLIC** Mölnlycke Health Care s.r.o.  
Hajkova 2747/22, 130 00 Prague 3  
Tel: +420 221 890 517 (**Reception**)  
Tel: +420 221 890 511 (**Customer Service**)

**DENMARK** Mölnlycke Health Care ApS,  
Gydevang 33, DK-3450 Allerød  
Tel: +45 48 16 82 68  
**Customer Service** Tel: +80 886 810

**ESTONIA** Mölnlycke Health Care OÜ,  
Vabaõhumuuseumitee 4-67, EE-13522 Tallin  
Tel: +372 671 1520

**FINLAND** Mölnlycke Health Care Oy,  
Pitäjänmäentie 14, 3 krs, FI-00380 Helsinki  
Tel: +358 201 622 300  
**Customer Service** Tel: 0800 113 442

**FRANCE** Mölnlycke Health Care SAS,  
13 Allée du Château Blanc  
59290 WASQUEHAL  
Tel: +33 320 122 555  
**Customer Service** Tel: 0800 910 292

**GERMANY** Mölnlycke Health Care GmbH  
Grafenberger Allee 297  
DE-40237 Düsseldorf  
Tel: +49 211 92 08 80  
**Customer Service** Tel: 0800 186 21 80

**HUNGARY** Mölnlycke Health Care Kft  
HU-1134 Budapest, Dévai u. 26-28. III. em.  
Tel: +36 1 477 3080

**INDIA** Mölnlycke Health Care India Pvt Ltd  
Universal Business Centre  
Second Floor, Universal Trade Tower  
Sector 49, Gurgaon Sohna Road  
Gurgaon 122 018, Haryana  
Tel: +91 124 4696 008

**ITALY** Mölnlycke Health Care s.r.l.,  
Via Marsala 40/C, IT-21013 Gallarate (VA)  
Tel: +39 0331 714 011

**JAPAN** Mölnlycke Health Care KK,  
4F Concieria Nishi-Shinjuku Tower's West  
6-20-7 Nishi-Shinjuku,  
Shinjuku-Ku, Tokyo 160-0023  
Tel: +81 3 6914 5004

**KOREA** Mölnlycke Health Care Korea Co., Ltd  
2nd floor, Hansol Richville.  
46 Wiryeseong-daero,  
Songpa-gu, Seoul, 05627  
South Korea  
Tel: +82 2 3789 1402

**LATVIA** Mölnlycke Health Care AB,  
Dzelzavas 120 M, Riga LV-1021  
Tel: +372 6711 522

**LITHUANIA** Mölnlycke Health Care AB,  
J. Galvydzio g. 3, LT 08236 Vilnius  
Tel: +370 5 2274 58 37

**NETHERLANDS** Mölnlycke Health Care B.V.  
Postbus 3196, 4800 DD Breda  
Tel: +31 76 521 9663

**NORWAY** Mölnlycke Health Care AS  
Postboks 6229 Etterstad, NO-0603 Oslo  
Tel: +47 22 70 63 70  
**Customer Service** Tel: 800 161 78

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Mölnlycke Health Care Ltd.
**USA** Mölnlycke Health Care US, LLC