

Stays on Uniquely conforms



Mepilex[®] Border Flex
Mepilex[®] Border Flex Lite




Mölnlycke[®]

Wounds need to heal. Patients need care. Nurses need time. Budgets are tight. The challenge is finding a dressing that is comfortable and stays on as long as it is needed.

Introducing the Mepilex® Border Flex range

Our next generation of flexible dressings are designed to stay on and uniquely conform aiming to give time back to nurses, cost savings back to managers and quality of life back to patients.




Track the spread of exudate without lifting dressing¹



Flex Technology™ to conform⁴



Trap exudate containing bacteria away from wound²



Channel exudate to maximise capacity⁵



Absorb exudate with minimal swelling³




Safetac® for gentle dressing changes⁶




*Not available in Mepilex Border Flex Lite

Choose the right dressing for your patient's needs


Designed for a wide range of exuding wounds including: leg and foot ulcers, traumatic wounds (e.g. skin tears) and surgical wounds, as well as part of a prophylactic therapy to help prevent skin damage, e.g. pressure ulcers.



Square




Oval



Rectangular

Mepilex® Border Flex Lite

Mepilex® Border Flex Lite, a four layer dressing with Flex Technology™ for non- to moderately exuding wounds also available in smaller sizes.



Square



Rectangular

More than half (58%) of patients complain that their dressing pulls or doesn't move with their body³.

Designed to conform and stay on


Thanks to Flex Technology™, the unique conformability of Mepilex® Border Flex allows it to adapt to the shape and movement of the patient¹³. Reducing pressure on the skin and pull on the borders increases comfort and minimises the risk of detachment⁶. And it's water resistant so your patients can shower⁷.



Unique Flex Technology™ = 360° flexibility



Y-cuts expand to allow movement:



Mepilex® Border Flex is as much as up to **11x** more flexible than other bordered foam dressings tested¹³.

Up to 1/3 of patients complain that their dressing needs changing daily³.


Smart exudate management

Giving you confidence to leave the dressing on for longer and still maintain an optimal wound healing environment.




1. Rapid absorption with minimal swelling:

The foam layer allows rapid absorption of both normal and viscous exudate^{7,8}.




2. Maximises capacity by channelling exudate:

Exudate is channelled away from the wound and over a wider area to increase moisture vapour loss and reduce swelling^{3,5}.



3. Traps bacteria away from the wound:

The superabsorbent fibres in the retention layer prevent exudate containing bacteria from re-entering the wound bed, even under compression⁹. Backing film prevents entry of external bacteria¹⁰.



4. Tracks exudate without dressing removal:

The Exudate Progress Monitor allows the spread of exudate to be objectively assessed and recorded, without disturbing the wound¹.


*Not available in Mepilex Border Flex Lite. **particles/ microbes larger than 25 micro meters

How the Exudate Progress Monitor works

The Exudate Progress Monitor gives you the confidence to leave the dressing in place for longer¹.




The Exudate Progress Monitor is a grid of equidistant dots.




e.g. 4x4cm

As exudate spreads you can record a dot count that reflects exudate volume in the dressing.



e.g. 8x10cm

You can then refer to the recorded dot count to help monitor exudate trends.



e.g. 10x14cm

Consider dressing change when exudate is approaching the borders of the dressing*.

*Wound inspection and dressing change frequencies are driven by clinical decision and should be at the discretion of the treating clinician

9/10 patients feel their wound impacts their quality of life³.


Reduced pain and damage to the skin

All Mepilex® Border Flex dressings seek to reduce trauma caused by dressing changes. This is made possible by our proven Safetac® soft silicone wound contact layer, which minimises pain and damage to the skin and risk of maceration compared with traditional (adhesive) dressings^{5,11,12}.

More than

9/10

patients preferred dressings with Safetac® over other dressings with traditional adhesives⁵



Skin stripping occurs with traditional adhesive



Safetac® minimises skin stripping

Mepilex® Border Flex can make a positive impact on your budget

Performed under real world conditions, two individual QIPs conducted in hospital setting in the US showed substantial reductions in dressing utilization compared to their previous formulary dressings* whilst maintaining the healing trajectory of the wounds evaluated. A change in protocol due to the extended wear time of Mepilex Border Flex thus resulted in fewer dressing changes, documented cost reductions and less nursing time^{14,15}.

Visit molnlycke.com to find out how Mepilex® Border Flex dressings can help your patients, your team and your budgets.

Mepilex Border Flex

Order number	Size (cm)	Size (inch)	Wound pad size (cm/inch)	Dressings per pack
595200	7.5 x 7.5	3 x 3	4.5 x 4.5/1.8 x 1.8	5
595300	10 x 10	4 x 4	6.5 x 6.5/2.6 x 2.6	5
595000	12.5 x 12.5	5 x 5	8.5 x 8.5/3.3 x 3.3	5
595400	15 x 15	6 x 6	11 x 11/4.3 x 4.3	5
595600	15 x 20	6 x 8	11 x 16/4.3 x 6.3	5

Note: Please contact your local Mölnlycke representative for information about order numbers available in your country.




Mepilex Border Flex Oval

Order number	Size (cm)	Size (inch)	Wound pad size (sq cm/sq inch)	Dressings per pack
583500	7.8 x 10	3.1 x 4	15/2.3	5
583300	13 x 16	5.1 x 6.3	85/13.2	5
583400	15 x 19	6 x 7.5	134/20.7	5



Mepilex Border Flex Lite

Order number	Size (cm)	Size (inch)	Wound pad size (cm/inch)	Dressings per pack
581011	4 x 5	1.6 x 2	2 x 3/0.8 x 1.2	10
581100	5 x 12.5	2 x 5	2.5 x 8.5/1 x 3.3	5
581200	7.5 x 7.5	3 x 3	4.5 x 4.5/1.8 x 1.8	5
581300	10 x 10	4 x 4	6.5 x 6.5/2.6 x 2.6	5
581500	15 x 15	6 x 6	11 x 11/4.3 x 4.3	5



*Allevyn® Gentle Border and Optifoam® Gentle Border SA

References: 1. Mölnlycke Health Care. Data on File. 2019-21. 2. Mölnlycke Health Care. Data on File-24. 3. Mölnlycke Health Care. Data on File. 2018-17. 4. Mölnlycke Health Care. Data on File. 2018-27. 5. White, R. A multinational survey of the assessment of pain when removing dressings. Wounds UK 2008;4(1):14-22. 6. Mölnlycke Health Care. Data on File. 2018-13. 7. Mölnlycke Health Care. Data on File-22. 8. Mölnlycke Health Care. Data on File. 2019-20. 9. Mölnlycke Health Care. Data on File-24. 10. Mölnlycke Health Care. Data on File. 2018-25. 11. Woo, K., Coutts, P.M., Price, P., Harding, K., Sibbald, R.G. A randomised crossover investigation of pain at dressing change comparing foam dressings. Advances in Skin and Wound Care 2009;22(7):304-310. 12. Meaume S, et al. A study to compare a new self-adherent soft silicone dressing with a self-adherent polymer dressing in stage II pressure ulcers. Ostomy Wound Management, 2003. 13. 2. Mölnlycke Health Care. Data on File. 2019-12. 14. Nelson, D. A New Bordered Foam Dressing Technology Improves Wound Outcomes and Satisfaction While Reducing Dressing Utilization in Acute Care. Poster Presentation at WOCN Annual Meeting, Nashville, 2019. 15. Tyson, L. Study First: Driving the Case for Improving Hospital Wound Care. Poster Presentation at WOCN Annual Meeting, Nashville, 2019.

Numbers after hyphen are connected to references in internal document

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, P.O. Box 6, Entreprenörstråket 21, SE-431 21 Mölndal Sweden. Phone + 46 31 722 30 00. The Mölnlycke, Mepilex and Safetac trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care group of companies. Flex Technology is a trademark of Mölnlycke Health Care. ©2025 Mölnlycke Health Care AB. All rights reserved. HQIM007251