Mölnlycke®

Granulox®

Haemoglobin Spray
Oxygenates wounds

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
Diapharm Switzerland GmbH
Hochbergerstrasse 70
CH-4057 Basel

CH Importer:
Mölnlycke Health Care AG
Brandstrasse 24
CH-8952 Schlieren

Australian sponsor address:
Mölnlycke Health Care Pty. Ltd.
12 Narabang Way, Belrose
NSW 2085 AUSTRALIA

Toll free number:
USA 1-800-882-4582
Canada 1-800-494-5134

en Medical Device	no Medisinsk utstyr
de Medizinprodukt	sl Medicinski pripomoček
fr Dispositif médical	bg Медицинско изделие
es Producto sanitario	ro Dispozitiv medical
nl Medisch hulpmiddel	sk Zdravotnícka pomôcka
sv Medicinteknisk produkt	tr Tıbbi cihaz
it Dispositivo medico	lt Medicinos priemonė
fi Lääkinnällinen laite	lv Medicīniska ierīce
pt Dispositivo médico	et Meditsiiniseade
da Medicinsk udstyr	ru Медицинское изделие
el ιατροτεχνολογικό προϊόν	hr Medicinski proizvod
pl Wyrób medyczny	sr Medicinsko sredstvo
cs Zdravotnický prostředek	is Lækningatæki
hu Orvostechnikai eszköz	ar جهاز طبي

Granulox®

Haemoglobin Spray. Oxygenates wounds.

Product description

Granulox® Haemoglobin Spray is a red, aqueous solution for treatment of chronic hypoxic wounds. Chronic wounds are usually characterized by a below-average oxygen supply to the affected tissue. Granulox provides the wound with the required oxygen by means of diffusion. The active substance, Haemoglobin, topically supplies the base of the wound with oxygen. The improved oxygen supply to the base of the wound supports wound healing.

Content: 10% carbonylated Haemoglobin, 0.7% phenoxyethanol, 0.9% sodium chloride, 0.05% N-acetylcysteine and water.

Indications for use

Granulox is used for treatment of chronic wounds, such as:

- venous leg ulcer
- arterial leg ulcer,
- mixed leg ulcer,
- diabetic foot ulcers,
- secondary healing of surgical wounds,
- pressure ulcers.

Granulox can also be used on sloughy and infected wounds after appropriate treatment.

Contraindications

Granulox should not be used on pregnant women since use in this condition has not been assessed. Granulox should not be used simultaneously with locally acting medicines, such as antibiotics, because no data on potential interactions with Haemoglobin is available.

Precautions

- Do not use Granulox on patient and/or user with known hypersensitivity to the product materials/ components.
- Some wound antiseptics are known to impair the efficacy of Granulox. Therefore, after application of an antiseptic, the wound area must be flushed thoroughly with a physiological solution.
- Do not use with occlusive plastic films and hydrocolloids.
- Side effects of Granulox have not been reported. If your symptoms worsen or you develop side effects, stop using Granulox and consult a doctor.

Instructions for use

Granulox can be used by health care professionals or lay persons under supervision of health care professionals and is suitable for repeated use.

The precondition for use is appropriate medical treatment of the primary disease. In the case of chronic arterial wounds, alternatives within best standard of care for improving blood flow must already have been exhausted.

Before starting treatment using Granulox, appropriate debridement and cleaning of the wound is required. Some wound antiseptics are known to impair the efficacy of Granulox. Therefore, after application of an antiseptic and after use of proteolytic (enzymatic) debridement, the wound area must be flushed thoroughly with a physiological solution.

Granulox should be sprayed onto the cleaned wound as follows:

- Remove the protective cap from the spray container.
- On first use, nonstop spraying may occur due to the very tight construction of the nozzle. If so, simply jiggle the nozzle gently back and forth and guide it slightly upwards. Spraying will then stop. When using Granulox for the first time, you may choose to spray into a tissue or over a sink.
- Hold the nozzle about 5 to 10 cm from the wound and spray until the wound is completely coated with a thin film. A spraying time of 1–2 seconds is normally sufficient to cover a wound area of 2×3 cm.
- After use, put the protective cap back on the spray container.
- If the nozzle becomes blocked after several uses, the spray head can be replaced with a new one supplied in the package. Alternatively, the spray head can be cleaned by using a sterile or preserved solution for cleaning. After cleaning, the first spray should be discarded.
- Apply the product every time the dressing is changed, or at least every three days.
- After use, cover the wound with a new permeable wound dressing.
- Granulox remains at a constant consistency, so there is no need to shake. And unlike conventional spray cans, shaking gives no indication of the amount of Granulox still remaining.
- One pack of Granulox is sufficient for up to 30 applications depending on the size of the wound. It empties 99% of contents.

After opening, Granulox can be used until the expiry date.

Special storage conditions and handling conditions

Granulox should be stored in a refrigerator at temperatures from 2°C to 8°C.

On the days it is being used, Granulox can be stored at room temperature (max. 25°C).

If Granulox is used on a daily basis, the spray container can be continuously stored at room temperature until it is empty. The storage time at room temperature must not exceed a maximum period of 6 weeks.

Keep out of reach of children. Pressurized container: May burst if heated. Do not pierce or burn, even after use.

Disposal should be handled according to local environmental procedures.

After use, dispose of the products as clinical waste in accordance with local regulations.

Other information

Any serious incident that occurs from the use of Granulox should be reported to the manufacturer and the competent authority in your country.

For EU-customers, further information on the safety and clinical performance of Granulox® with Basic UDI-DI: 42517781MD910DE, please refer to the 'Summary of Safety and Clinical Performance for Granulox®'. You can access the document on <https://ec.europa.eu/tools/eudamed>.

Granulox® is registered trademark of Mölnlycke Health Care AB.

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