

Mölnlycke® Z-Flo™ Fluidized positioner – continuous protection

Adult and paediatric



Patient positioning tools are inadequate

Patient positioning is a cornerstone of pressure ulcer prevention and treatment. And yet the tools that you have at your disposal is less than adequate.

Ask yourself if blankets, pillows and wedges offers the following benefits:

- Pressure redistribution through envelopment?
- Ability for you to dictate the desired position?
- Ability to hold that position over time?
- Ability to fold it away to protect sensitive areas?

Continuous protection by pressure redistribution and lasting positioning

Mölnlycke® Z-Flo™ Fluidized positioner allows you to mold it to fit the patient's anatomy for pressure redistribution. It allows you to tailor it to individual positioning needs and will maintain the desired position over time.



Areas in the positioner **can be depressed, or hollowed out**, for medical lines, tubes or existing wounds.



Place Z-Flo under the head and align it with the top of the shoulders



Mold Z-Flo so that it conforms to the shape of the head and fills in the cervical spine area



Positioner in place offering both **positioning and pressure relief**

Mölnlycke® Z-Flo™ in the ICU

Mölnlycke Z-Flo is a versatile product with multiple uses. In the ICU, the of the most common ones is for occipital pressure injury prevention and patient turning.

Reduction by
87.7%
in occipital
pressure injuries!

Occipital pressure injuries prevention

One study sought to evaluate the clinical efficacy and feasibility of the Mölnlycke Z-Flo Fluidized positioner in reducing the occurrence of occipital PIs, in an ICU setting. They choose to study patients who were receiving ECMO where frequent off-loading and repositioning of the occiput is a challenge¹.

Utilize Mölnlycke Z-Flo to increase the contact surface area and for repositotining of the head to aid in your pressure injury prevention program.




Turning and repositioning

Mölnlycke Z-Flo Fluidized positioner will maintain its shape – and the patient's position over time.

In addition, the fluidized positioner lets the clinician dictate the degree of turn, ranging from a small microturn all the way up to large-degree turns.

One study sought to compare using Mölnlycke Z-Flo Fluidized positioner with using pillows. Results showed a statistically significant difference between the two positioning devices' ability to maintain a turn over time².



Proving it every day

At Mölnlycke, we deliver innovative solutions for managing wounds, improving surgical safety and efficiency and preventing pressure ulcers. Solutions that help achieve better outcomes and are backed by clinical and health-economic evidence.

In everything we do, we are guided by a single purpose: to help healthcare professionals perform at their best. And we're committed to proving it every day.

Adult Fluidized positioner

Product Code	Description	Dimensions	Psc/Case
1401007	Small positioner	29 x 50 cm	8
1401011	Medium positioner	40 x 76 cm	1
1401005	Large positioner	64 x 91 cm	1

Paediatric Fluidized positioner

Product Code	Description	Dimensions	Psc/Case
1400239-10	Medium positioner	23 x 38 cm	6
1400200-10	Large positioner	30 x 51 cm	6



Learn more how Z-Flo™ Fluidized positioner can help you and your patients

References: 1. Barakat-Johnson M et al. Evaluation of a fluidised positioner to reduce occipital pressure injuries in intensive care patients: A pilot study. *Int Wound J.* 2018;1-9.
2. Kapp S, Gertz M, Gefen A, Prematunga R, Santamaria N. An observational study of the maintenance of the 30°side-lying lateral tilt position among aged care residents at risk of developing pressure injuries when using the standard care pillow and a purpose-designed positioning device. *Int Wound J.* 2019.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00.
The Mölnlycke trademark, name and logo are registered globally to one or more of the Mölnlycke Health Care group of companies. Z-Flo is a trademark of Mölnlycke Health Care ©2022 Mölnlycke Health Care AB. All rights reserved. HQIM003597

