Mepilex® Lite

Instruction for Use

Product description

Mepilex Lite is an absorbent self-adhesive soft silicone foam dressing, its composition is as follows:

A layer of soft silicone wound contact layer, a layer of polyurethane foam absorption layer, a layer of breathable and waterproof outer membrane.

Safetac® technology

Safetac® it is a technology of soft silicone, which can reduce the patient's pain and trauma to the wound. The reason why Safetac® technology is less painful is because it:

1. Gently stick on a dry surface, such as skin, but will not stick on a wet surface, such as an open wound;
2. Covers skin pores and more skin surface, thereby preventing it from being peeled off on the skin;
3. Sealing the wound area, to reduce risk of maceration.

Mode of action

Mepilex Lite is a highly conformable dressing that absorbs exudates, maintains a moist wound environment.

The Safetac layer will seal the edges of the wound, preventing exudates from leaking into the surrounding normal skin, reducing the risk of maceration.

Mepilex Lite can be cut according to the shape and position of the wound.

Mepilex Lite can be combined with the use of bandages.

Intended use

Mepilex Lite can be used for the management of non to low exuding wounds such as leg and foot ulcers, pressure ulcers, partial thickness burns, radiation skin reactions and Epidermolysis Bullosa. Mepilex Lite can also be used as a protection of compromised and/or fragile skin.

Instruction for use

1. Cleanse the wound in accordance with clinical practice.
2. Thoroughly dry the skin around the wound.
3. Select an appropriate dressing size. The dressing should cover the dry surrounding skin by at least 1–2 cm for small sizes (sizes up to 12,5×12,5 cm) and 3–5 cm for large sizes. If required, the dressing may be cut to suit various wound shapes and locations.
4. Remove the first release film and apply the adherent side to the wound.
5. Remove the remaining release film and smooth down the dressing on the skin.
Do not stretch the dressing.

6. When necessary, fixate the dressing with a bandage or other fixation.

Replacement frequency

Mepilex Lite can be kept on the wound for several days according to the wound condition and surrounding skin condition, or as indicated by clinical practice.

A change in dressing regimen may result in an initial increased level of exudate that temporarily may require an increased frequency of dressing changes.

Disposal should be handled according to local environmental procedures.

Precautions

- Do not use on patients with known hypersensitivity to the dressing or its materials.
- Do not use together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- If you see signs of clinical infection e.g. fever or if the wound or surrounding skin becomes red, warm or swollen, consult a health care professional for appropriate treatment.
- Mepilex Lite is ethylene oxide sterilization of aseptic packaging, such as packaging is damaged or opened prior to use, do not use. This product cannot be sterilized again for use.
- Do not reuse, if reused performance of the product may deteriorate and cross contamination may occur.

Storage

Mepilex Lite should be stored in a dry environment below 35°C (95°F) and protected from light.

Sterilization

Sterile packaging, sterilized with ethylene oxide. Do not use if sterile barrier is damaged or opened prior to use.

Do not re-sterilise.

Do not use expired products.

Shelf-life: Three years.

Other information

If any serious incident has occurred in relation to the use of Mepilex Lite, it should be reported to Mölnlycke Health Care and to your local competent authority.

Mepilex® and Safetac® are registered trademarks of Mölnlycke Health Care AB.

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Telephone: +46 31 722 30 00

Fax: +46 31 722 34 00

Agent and after-sales service agency contact information

Name: Mölnlycke Health Care (Shanghai) Co., Ltd

Address: Room 629, No 8, Huajing road, China (Shanghai) pilot free trade zone

Zip code: 200131

Telephone: 010-51288571; 021-62366898

Fax: 010-84554552; 021-62366830

Website: www.molnlycke.cn

Registration License Number: 国械注进 20173141524

Product Technical Requirements No.: 国械注进 20173141524

Chinese IFU version date: 11th Jun, 2025

Chinese IFU version: 20171524 REV.6WF

See the Chinese label for other content

Symptom indication:

	Do not re-use		Catalogue number
	Keep dry		Unique Device Identifier
	Expiry date		Sterilize Sterilization method: ethylene oxide sterilization
	Batch code		Single sterile barrier system
	Manufacturer		Do not use if package is damaged
	Medical device		Caution. See instructions for use
	CE-mark with Notified Body identification number Compliant with European Medical Device Regulation (MDR)		

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