Mepitel® One Soft Silicone Wound Contact Layer

Instructions for use

Product Description

Mepitel One is a transparent and soft wound contact dressing, with one side composed of the porous polyamide film applied with adherent silicone (Safetac).

Safetac® technology

Safetac® is a soft silicone adhesive technology that reduces pain and trauma to the wound surface. Safetac® technology is less painful because it:

1. Can be gently adhered to dry surfaces such as skin, but not to wet surfaces such as open wounds;
2. Covers skin pores and more skin surface to prevent peeling from the skin;
3. Sealing the wound area to ensure that exudate does not spread to the surrounding skin and reduce the risk of maceration.

Mode of Action

- Mepitel One is not absorbent. The open perforated structure allows exudate to pass vertically into a secondary absorbent dressing pad, which must be changed as required by the conditions of the wound and the amount of exudate in order to prevent maceration.
- The integrity of the Mepitel One dressing reduce the necessity for frequent primary dressing changes and allows secondary dressing changes with minimised pain as required. Due to the non-sticky top surface and adequate adhesion to skin it is also possible to use Mepitel One as a stand alone product for protection of damaged skin.
- Mepitel One can be used under compression bandages.
- Mepitel One can be cut to different sizes to suit the shape and location of the wound.

Intended use

Mepitel One is a wound contact layer dressing designed for the treatment of exudative wounds on the body surface such as surgical incisions, partial thickness burns, diabetic foot and surgical trauma (abrasions, lacerations).

Instructions for use

Mepitel One with this symbol can help you apply and remove the dressing on skin tears. For these wounds, there is a risk of re-opening the flap when removing the dressing. Apply the dressing with the arrow pointing in the

same direction as the flap. Start application from the base of the flap. When removing the dressing, begin removal following the printed symbol, in the direction the arrow is pointing.

1. Cleanse the wound in accordance with clinical practice and dry the surrounding skin thoroughly.
2. Choose a size of Mepitel One that covers the wound and the surrounding skin by at least 2 cm. For larger wounds, more overlap is required. If more than one piece of Mepitel One is used, overlap the dressings. Make sure that the holes are not blocked.
3. Remove the protective plastic layers and apply Mepitel One with the sticky side to the wound.
4. The dressing is applied in a correct way when you can read the text printed on the dressing. Smooth the dressing in place onto the surrounding skin to get a good seal.
5. Apply an outer absorbent dressing, on top of Mepitel One and fixate.

Dressing change

Mepitel One may be left in place for up to 7 days depending on the condition of the wound and surrounding skin, or as indicated by clinical practice. To prevent maceration, exudate should pass freely through the dressing and the holes should not be blocked.

If the absorption of the secondary dressing pad reaches saturation, it should be replaced, while the Mepitel One does not need to be replaced.

The Mepitel One can be used independently.

Warning

Mepitel One should be used correctly to avoid leaving marks when used for mesh skin grafting after burns or facial plastic surgery.

Precautions

- If you see signs of infection e.g. fever or the wound or surrounding skin becoming red, warm or swollen, consult a health care professional for appropriate treatment.
- If Mepitel One is used on Epidermolysis Bullosa patients, employ extra surveillance at dressing changes. When used on partial thickness burns with high risk of rapid granulation, avoid placing pressure on the dressing.
- When used on after facial resurfacing, pressure should be avoided when removing the dressing, lift and reposition the dressing at least every second day.
- When used on bleeding wounds or wounds with high viscosity exudate, Mepitel

One should be covered with a moist absorbent dressing pad.

- When using Mepitel One for the fixation of skin grafts and protection of blisters, the dressing should not be changed before the fifth day post application.
- Determine the replacement frequency based on physician's clinical procedures.

Storage

Mepitel One should be stored in a dry environment at room temperature.

Sterilization

Mepitel One is sterilized by ethylene oxide. Do not use if inner package is damaged or opened prior to use.

Do not re-sterilize.

Product model/ specification

Product model	Specification
289000	5×7cm
289100	5×7.5cm
289170	6×7cm
289200	8×10cm
289270	9×10cm
289300	7.5×10cm
289400	12×15cm
289470	13×15cm
289500	10×18cm
289700	17×25cm
289670	24×27.5cm

Period of validity: 3 years.

Manufacturer

Manufacturer: Mölnlycke Health Care AB

Manufacture's Address: Entreprenörsstråket 21 SE-431 53 Mölndal Sweden

Address of Manufacturer Site: Postfach 76 SF-50101 Mikkeli 10 Finland

Made in Finland

Tel: +46 31 722 30 00

Fax: +46 31 722 34 00

Legal agent/after-sale service agent in China mainland

Legal agent: Mölnlycke Health Care (Shanghai) Co., Ltd.

Address: Room 629, No. 8, Huajing Road, China (Shanghai) Pilot Free Trade Zone

Postcode: 200131

Tel: 010-51288571; 021-62366898

Revision No.: 20161486 REV.7WF

Fax: 010-84554552; 021-62366830

Website: www.molnlycke.cn

Registration number: 国械注进20162141486

Technical Requirement No.: 国械注进20162141486

Production date/expiration date: see label.

Revision date: 3th Jun 2025

Revision No.: 20161486 REV.7WF

Logo description:

	Batch code		Do not re-use
	Expiry date		Sterilize Sterilization method: ethylene oxide sterilization
	Catalogue number		Single sterile barrier system
	Manufacturer		Do not use if package is damaged
	Medical device		Caution. See instructions for use
	CE-mark with Notified Body identification number Compliant with European Medical Device Regulation (MDR)		

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