

# EU Type Examination Certificate

This is to certify that:

Mölnlycke Health Care AB  
Gamlestadsvägen 3C  
Box 13080  
Göteborg  
SE-402 52  
Sweden

Holds Certificate Number:

CE 687999

In respect of:

**Biogel® gloves for Personal Protection.**  
**Various models, see continuation sheets for details.**

on the basis that BSI carried out the relevant Type Examination procedures under the requirements with the Regulation (EU) 2016/425 of the European Parliament and Council relating to Personal Protective Equipment Regulation (PPE) Annex V (Module B) and meets the relevant health and safety requirements specified in Annex II

For and on behalf of BSI, a Notified Body for the above Regulation (Notified Body Number 2797):

Previous Notified Body: BSI 0086

First Issued: 2018-08-21

Latest Issue: 2019-09-12



Drs. Dave Hagenaaars, Managing Director

Effective Date: 2019-09-12

Expiry Date: 2023-08-21

Page: 1 of 9



...making excellence a habit.™

# EU Type Examination Certificate

No. CE 687999

## Product Specification

**PPE Category:** Category III - Complex.

**Applicable Standards:** The following Harmonized European Standards:

EN 420:2003+A1:2009 - Protective gloves. General requirements.

EN 374-1:2016 - Protective gloves against chemicals and micro-organisms - Terminology and performance requirements.

EN 374-2:2014 - Protective gloves against dangerous chemicals and micro-organisms - Determination of resistance to penetration.

EN 374-4:2013 – Protective gloves against chemicals and micro-organisms. Determination of resistance to degradation by chemicals.

EN 374-5:2016 – Protective gloves against dangerous chemicals and micro-organisms. Terminology and performance requirements for micro-organism risks.

EN 16523-1:2015 - Protective gloves against dangerous chemicals and micro-organisms. Determination of resistance to permeation by chemicals.

EN 16604:2004 – Clothing for protection against contact with blood and body fluids. Determination of resistance of protective clothing materials to penetration by blood-borne pathogens. Test method using Phi-X174 Bacteriophage.

First Issued: 2018-08-21

Latest Issue: 2019-09-12

Effective Date: 2019-09-12

Expiry Date: 2023-08-21

Page: 2 of 9

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.

# EU Type Examination Certificate

No. CE 687999

**Model:** Biogel® NeoTech S44509 (formerly Skinsense®)

**Product Code:** S44509XX  
Size range 5.5 to 9

**Classification:** Protective gloves for use against micro-organisms and general applications. Designed for single use by individuals including those who are sensitive to natural rubber latex to provide cover to the hands and wrists.

**Description :** A sterile, powder-free, synthetic polychloroprene surgical overglove with polymeric inner surface treatment to assist donning. The glove provides a barrier to protect against cross-infection.

## EN 420:2003+A1:2009 - General requirements:

Characteristic	Performance
pH value	Pass
Dexterity	Level 5
Length and Fit	Pass

## EN 374-2:2014 – Resistance to penetration

Performance Level 3

## EN 374-4:2013 – Resistance to degradation

Tested for all chemicals listed below.

**EN 374-1:2016 Resistance to chemical permeation to test method**  
**EN 16523-1:2015 Determination of resistance of permeation by chemicals.**  
Type C

Chemical	Permeation Level
40% Sodium Hydroxide (K)	6

## EN 16604:2004 Resistance to penetration by blood-borne pathogens

Pass

First Issued: 2018-08-21  
Latest Issue: 2019-09-12

Effective Date: 2019-09-12  
Expiry Date: 2023-08-21

Page: 3 of 9

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.

# EU Type Examination Certificate

No. CE 687999

**Model:** Biogel® Tech  
**Product Code:** S44822XX  
**Size range 5.5 to 9**

**Classification:** Protective gloves for use against micro-organisms and general applications. Designed for single use and to provide cover to the hands and wrists, provides low chemical protection.

**Description :** A sterile, powder-free, natural rubber latex surgical glove with polymeric inner surface treatment. The glove is anatomically shaped, has curved fingers and a smooth finish.

## EN 420:2003+A1:2009 - General requirements:

Characteristic	Performance
pH value	Pass
Dexterity	Level 5
Length and Fit	Pass
Protein Content	<50µg/g

## EN 374-2:2014 – Resistance to penetration

Performance Level 3

## EN 374-4:2013 – Resistance to degradation

Tested for all chemicals listed below.

**EN 374-1:2016 Resistance to chemical permeation to test method**  
**EN 16523-1:2015 Determination of resistance of permeation by chemicals.**  
**Type C**

Chemical	Permeation Level
40% Sodium Hydroxide (K)	6

## EN 16604:2004 Resistance to penetration by blood-borne pathogens

Pass

First Issued: 2018-08-21  
Latest Issue: 2019-09-12

Effective Date: 2019-09-12  
Expiry Date: 2023-08-21

Page: 4 of 9

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.

# EU Type Examination Certificate

No. CE 687999

**Model:** Biogel® Tech  
**Product Code:** S44961XX  
**Size range 5.5 to 9**

**Classification:** Protective gloves for use against micro-organisms and general applications. Designed for single use and to provide cover to the hands and wrists.

**Description :** A sterile, powder-free, natural rubber latex surgical glove with polymeric inner surface treatment. The glove is anatomically shaped, has straight fingers and a micro rough finish.

## EN 420:2003+A1:2009 - General requirements:

Characteristic	Performance
pH value	Pass
Dexterity	Level 5
Length and Fit	Pass
Protein Content	<50µg/g

## EN 374-2:2014 – Resistance to penetration

Performance Level 3

## EN 374-4:2013 – Resistance to degradation

Tested for all chemicals listed below.

**EN 374-1:2016 Resistance to chemical permeation to test method**  
**EN 16523-1:2015 Determination of resistance of permeation by chemicals.**  
**Type C**

Chemical	Permeation Level
40% Sodium Hydroxide (K)	6

## EN 16604:2004 Resistance to penetration by blood-borne pathogens

Pass

First Issued: 2018-08-21  
Latest Issue: 2019-09-12

Effective Date: 2019-09-12  
Expiry Date: 2023-08-21

Page: 5 of 9

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.

# EU Type Examination Certificate

No. CE 687999

**Model:** Biogel® PI Tech S44409 (formerly PI Ultratouch)

**Product Code:** S44409XX  
Size range 5.5 to 9

**Classification:** Protective gloves for use against micro-organisms and general applications. Designed for single use by individuals including those who are sensitive to natural rubber latex to provide cover to the hands and wrists.

**Description :** A sterile, powder-free, transparent synthetic polyisoprene surgical overglove with polymeric inner surface treatment to assist donning. The glove provides a barrier to protect against cross-infection.

## EN 420:2003+A1:2009 - General requirements:

Characteristic	Performance
pH value	Pass
Dexterity	Level 5
Length and Fit	Pass

## EN 374-2:2014 – Resistance to penetration

Performance Level 3

## EN 374-4:2013 – Resistance to degradation

Tested for all chemicals listed below.

**EN 374-1:2016 Resistance to chemical permeation to test method**  
**EN 16523-1:2015 Determination of resistance of permeation by chemicals.**  
Type C

Chemical	Permeation Level
40% Sodium Hydroxide (K)	6

## EN 16604:2004 Resistance to penetration by blood-borne pathogens

Pass

First Issued: 2018-08-21

Latest Issue: 2019-09-12

Effective Date: 2019-09-12

Expiry Date: 2023-08-21

Page: 6 of 9

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.

# EU Type Examination Certificate

No. CE 687999

**Model:** Biogel® PI Tech Indicator Underglove

**Product Code:** S44416XX  
Size range 5.5 to 9

**Classification:** Protective gloves for use against micro-organisms and general applications. Designed for single use by individuals including those who are sensitive to natural rubber latex to provide cover to the hands and wrists.

**Description :** A sterile, powder-free, blue synthetic polyisoprene surgical underglove with polymeric inner surface treatment to assist donning. The glove provides a barrier to protect against cross-infection.

## EN 420:2003+A1:2009 - General requirements:

Characteristic	Performance
pH value	Pass
Dexterity	Level 5
Length and Fit	Pass

## EN 374-2:2014 – Resistance to penetration

Performance Level 3

## EN 374-4:2013 – Resistance to degradation

Tested for all chemicals listed below.

**EN 374-1:2016 Resistance to chemical permeation to test method**  
**EN 16523-1:2015 Determination of resistance of permeation by chemicals.**  
Type C

Chemical	Permeation Level
40% Sodium Hydroxide (K)	6

EN 16604:2004 Resistance to penetration by blood-borne pathogens

Pass

First Issued: 2018-08-21  
Latest Issue: 2019-09-12

Effective Date: 2019-09-12  
Expiry Date: 2023-08-21

Page: 7 of 9

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.

# EU Type Examination Certificate

No. CE 687999

## Certificate Administration Details

Technical File Reference: PPE Technical Files - Biogel® NRL Critical Environment Gloves and Biogel® Synthetic Critical Environment Gloves.

## Certificate Administration Record and BSI internal Technical File Review reference


Issue date	Comments	BSI Project Ref.
August 2018	First Issue	0086:18:8886142

## Certificate validity

The Certificate holder is responsible for ensuring that the Notified Body is advised of changes to any aspect of the overall process utilised in the manufacture of the products, failure to do so could invalidate the Certificate in respect of product manufactured following the introduction of such changes

The requirements for internal production control plus supervised product checks at random intervals, are referenced in BSI issued Certificate Number CE 578631.

The validity of the Certificate is also dependent on the maintenance of a documented quality system.



First Issued: 2018-08-21  
Latest Issue: 2019-09-12

Effective Date: 2019-09-12  
Expiry Date: 2023-08-21

Page: 8 of 9

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.



# EU Type Examination Certificate

No. CE 687999

## Location

Mölnlycke Health Care AB  
Gamlestadsvägen 3C  
Box 13080  
Göteborg  
SE-402 52  
Sweden

## Certified Activities

Biogel® gloves for Personal Protection.  
Various models, see continuation sheets for details.



First Issued: 2018-08-21  
Latest Issue: 2019-09-12

Effective Date: 2019-09-12  
Expiry Date: 2023-08-21

Page: 9 of 9

This certificate has been issued by and remains the property of BSI Group The Netherlands B.V., John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands and should be returned immediately upon request.  
To check its validity telephone +31 20 3460780. An electronic certificate can be authenticated [online](#).

BSI Group The Netherlands B.V., registered in the Netherlands under number 33264284, at John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands  
A member of BSI Group of Companies.