

# Tubigrip®

Multi-purpose tubular support bandage

Provides firm support for sprains, strains and swelling



## What size to choose?

Tubigrip size	Limb circumference (cm)	Limb circumference (inches)	Sub-bandage pressure (mmHg)
A	10 – 14.5	3.9 – 5.7	20 – 30
B	14.5 – 34	5.7 – 13.4	10 – 20
C	22 – 36	8.7 – 14.2	
D	26 – 41	10.2 – 16.1	
E	29 – 46	11.4 – 18.1	
F	33 – 53	13.0 – 20.9	
G	39 – 62	15.4 – 24.4	5 – 15
J	54 – 82	21.3 – 32.3	
K	70 – 113	27.6 – 44.5	
L	106 – 149	41.7 – 58.7	
M	119 – 167	46.9 – 65.7	

Tubigrip is available under a full range of sizes to suit most applications. The sub-bandage pressures are obtained in double layer.

Reference: Mölnlycke Health Care. Data on file.

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00. The Mölnlycke, Tubigrip trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care group of companies. © 2020 Mölnlycke Health Care AB. All rights reserved. HQIM002264.

