

Biogel® Synthetic Polyisoprene Glove Range



This product is not made from natural rubber latex.

Synthetic Polyisoprene

A sterile, powder-free, synthetic polyisoprene glove that eliminates the possibility of glove-related latex protein sensitisation because it is made from a synthetic elastomer and not natural rubber latex. The glove provides levels of fit, feel and comfort comparable to natural rubber latex because synthetic polyisoprene shares a similar molecular structure to natural rubber latex.



Biogel® PI Micro

Biogel® PI Micro Indicator® Underglove

Biogel® PI UltraTouch®

Biogel® PI Indicator® Underglove

Description	Enhanced tactile sensitivity	For use alone or as part of a Puncture Indication System	Enhanced sensitivity and softness	For use alone or as part of a Puncture Indication System
Color	Straw Color	Blue for Puncture Indication	Straw Color	Blue for Puncture Indication
Reorder Number	48555-48590	48955-48990	40955-40990	41655-41690
Thickness	8.3 mils*	8.3 mils*	10.6 mils*	11.0 mils*
Grip	1.5	1.0	1.0	1.0
Finish	Micro-roughened Surface	Smooth Finish for Double Donning	Micro-roughened Surface	Smooth Finish for Double Donning
Former	Curved Former	Curved Former	Curved Former	Curved Former
Cuff	Anti-Slip, Beaded Cuff	Beaded Cuff	Anti-Slip, Beaded Cuff	Beaded Cuff
Sizes	5.5-9.0	5.5-9.0	5.5-9.0	5.5-9.0
Quantities	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9
Tested for use with Chemotherapy agents	Yes	Yes	Yes	Yes
Ideal Use	For use in all surgeries where tactile sensitivity is desired.	For use in all surgeries. Specially designed for ease in double-gloving and enhanced tactile sensitivity.	For general purpose use in all surgeries.	For use in all surgeries. Specially designed for ease in double-gloving.

Biogel® Latex Glove Range



Natural Rubber Latex

A sterile, powder-free, natural rubber latex glove. All Biogel® latex gloves have exceptionally low extractable protein levels and low endotoxin levels (<0.5 EU/mL).



Biogel® Surgeons

Biogel® Eclipse®

Biogel® Eclipse® Indicator® Underglove

General purpose barrier protection	30% softer [‡] with enhanced sensitivity	For use alone or as part of a Puncture Indication System
Straw Color	Straw Color	Green for Puncture Indication
82255-82290	75155-75190	73255-73290
10.6 mils*	9.4 mils*	9.4 mils*
1.0	2.0	1.0
Micro-roughened Surface	Micro-roughened Surface	Smooth Finish for Double Donning
Curved Former	Curved Former	Curved Former
Beaded Cuff	Anti-Slip, Beaded Cuff	Beaded Cuff
5.5-9.0	5.5-9.0	5.5-9.0
50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9
Yes	Yes	Yes
For general purpose use in all surgeries.	For general purpose use in all surgeries especially when enhances tactile sensitivity is desired.	For use in all surgeries. Specially designed for ease in double-gloving.

Biogel® Synthetic Polychloroprene Glove Range



This product is not made from natural rubber latex.

Synthetic Polychloroprene

A sterile, powder-free, synthetic polychloroprene glove that reduces the possibility of glove-related latex protein sensitisation because it is made from a synthetic elastomer.



Biogel® NeoDerm®

Biogel® Skinsense®

Biogel® Skinsense® Indicator® Underglove

Made without chemical accelerators or CPC [†]	Provides extra protection due to longer chemical breakthrough time [‡]	For use alone or as part of a Puncture Indication System
Khaki Color	Straw Color	Blue for Puncture Indication
42955-42990	50955-50990	40655-40690
8.7 mils*	8.3 mils*	8.3 mils*
1.5	1.5	1.0
Micro-roughened Surface	Micro-roughened Surface	Smooth Finish for Double Donning
Curved Former	Curved Former	Curved Former
Beaded Cuff	Beaded Cuff	Beaded Cuff
5.5-9.0	5.5-9.0	5.5-9.0
50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9	50 Pairs per box/4 boxes per Case 40 pairs per box for size 9
No	Yes	Yes
For use in all surgeries especially for users needing to avoid chemical accelerators. [†]	For use in all surgeries or where extended protection against chemicals or cytotoxic drugs is required.**	For use in all surgeries or where extended protection against chemicals or cytotoxic drugs is required.** Specially designed for ease in double-gloving.

* Finger Thickness Single Wall, typically

** Breakthrough testing times are available by request

† Made from synthetic chloroprene excluding chemical accelerators known to cause skin irritation, such as Thiadiazoles, Thiurams, Carbamates, Thioureas and Diphenylguanidine.¹

‡ Compared to Biogel® Surgeons®