

# Latex Surgical Gloves

## Biogel Eclipse

A sterile, latex surgical glove for tactile sensitivity. It is 30% softer than our standard surgical gloves for improved fit, feel and comfort.



Ref	Size	Pairs	Ref	Size	Pairs
75155	5.5	50 / Box	75175	7.5	50 / Box
75160	6	50 / Box	75180	8	50 / Box
75165	6.5	50 / Box	75185	8.5	50 / Box
75170	7	50 / Box	75190	9	40 / Box

## Biogel Eclipse® Indicator® Underglove

A green, sterile, latex underglove that can be worn together with the Biogel Eclipse overglove to create a coloured Puncture Indication System.



Ref	Size	Pairs	Ref	Size	Pairs
73255	5.5	50 / Box	73275	7.5	50 / Box
73260	6	50 / Box	73280	8	50 / Box
73265	6.5	50 / Box	73285	8.5	50 / Box
73270	7	50 / Box	73290	9	40 / Box

## Biogel Eclipse® Indicator® System

A green latex indicator underglove and a straw-coloured overglove, creating a Puncture Indicator System proven to provide clear, fast and large coloured puncture indication<sup>1</sup>.



Ref	Size	Pairs	Ref	Size	Pairs
60755	5.5	2 x 25 /Box	60775	7.5	2 x 25 /Box
60760	6	2 x 25 /Box	60780	8	2 x 25 /Box
60765	6.5	2 x 25 /Box	60785	8.5	2 x 25 /Box
60770	7	2 x 25 /Box			

## Biogel® Surgeons

A sterile, latex surgical glove suitable for all procedures. The unique Biogel coating provides great fit, feel and comfort and makes the glove easy to don, even with damp hands.



Ref	Size	Pairs	Ref	Size	Pairs
82255	5.5	50 / Box	82275	7.5	50 / Box
82260	6	50 / Box	82280	8	50 / Box
82265	6.5	50 / Box	82285	8.5	50 / Box
82270	7	50 / Box	82290	9	40 / Box

Puncture indication. **Clear. Fast. Large.**



### Clear

The inner glove colour physics are based on science and selected to trigger sensitivity in the human eye. The outer gloves are opaque, so when punctured they alert you by showing the colour of the underglove<sup>1</sup>.

### Fast

The glove surface treatment highlights the puncture spot immediately by spreading trapped fluid to create a visible indication that grows quickly<sup>1</sup>.

### Large

The indication spot grows large so you can see it quickly - and so others can see it from a distance in case it's outside your direct line of sight<sup>1</sup>.

You can count on finding your preferred fit, feel and comfort.

<sup>1</sup>REPRO355, Spectrophotometer Measurement of the Colour Indicator Type Surgical Gloves.  
<sup>\*\*</sup>Breakthrough testing data available upon request REPR 0618

Find out more at:  
[www.molnlycke.com.au](http://www.molnlycke.com.au) or  
[www.molnlycke.co.nz](http://www.molnlycke.co.nz)


Molnlycke Health Care Pty Ltd, Level 4, 12 Narabang Way, Betrose, NSW 2085 T 02 9485 1544, F 02 9485 1555.  
 Molnlycke Health Care Pty Ltd, PO Box 62-027, Mt Wellington 1641, New Zealand. Orders & Enquiries 0800 005 231.

The Molnlycke, Biogel and Biogel Indicator trademarks, names and logos are registered globally to one or more of the Molnlycke Health Care group of companies. © 2026 Molnlycke Health Care AB. All rights reserved. AU\_002227



Biogel® - where your need for protection begins

Quick reference guide



## Synthetic Surgical Gloves

### Biogel® PI Micro

A sterile, synthetic surgical glove for extra tactile sensitivity that delivers the fit, feel and comfort of natural rubber latex. It is our thinnest surgical glove – it is 20% thinner than our regular Biogel PI gloves.



Ref	Size	Pairs	Ref	Size	Pairs
48555	5.5	50 / Box	48575	7.5	50 / Box
48560	6	50 / Box	48580	8	50 / Box
48565	6.5	50 / Box	48585	8.5	50 / Box
48570	7	50 / Box	48590	9	40 / Box

### Biogel® PI Micro Indicator® Underglove

A blue, sterile, synthetic underglove that can be used with a Biogel PI Micro overglove to create a coloured Puncture Indication System.



Ref	Size	Pairs	Ref	Size	Pairs
48955	5.5	50 / Box	48975	7.5	50 / Box
48960	6	50 / Box	48980	8	50 / Box
48965	6.5	50 / Box	48985	8.5	50 / Box
48970	7	50 / Box	48990	9	40 / Box

### Biogel® PI Micro Indicator® System

A sterile, double packed synthetic surgical glove that provides extra tactile sensitivity even when double gloving. It consists of a blue surgical underglove combined with a Biogel PI Micro overglove.



Ref	Size	Pairs	Ref	Size	Pairs
48355	5.5	2 x 25 /Box	48375	7.5	2 x 25 /Box
48360	6	2 x 25 /Box	48380	8	2 x 25 /Box
48365	6.5	2 x 25 /Box	48385	8.5	2 x 25 /Box
48370	7	2 x 25 /Box			

## Synthetic Surgical Gloves

### Biogel® PI UltraTouch®

A sterile, synthetic surgical glove suitable for all procedures that can be used alone or as part of the Biogel Puncture Indication System.



Ref	Size	Pairs	Ref	Size	Pairs
40955	5.5	50 / Box	40975	7.5	50 / Box
40960	6	50 / Box	40980	8	50 / Box
40965	6.5	50 / Box	40985	8.5	50 / Box
40970	7	50 / Box	40990	9	40 / Box

### Biogel PI Indicator® Underglove

A blue, sterile, synthetic underglove that can be used with a Biogel PI overglove to create a coloured Puncture Indication System.



Ref	Size	Pairs	Ref	Size	Pairs
41655	5.5	50 / Box	41675	7.5	50 / Box
41660	6	50 / Box	41680	8	50 / Box
41665	6.5	50 / Box	41685	8.5	50 / Box
41670	7	50 / Box	41690	9	40 / Box

### Biogel PI Indicator® System

A sterile, double packed synthetic surgical glove that consists of a blue underglove with a Biogel PI Overglove that can be worn together to create a coloured Puncture Indication System.



Ref	Size	Pairs	Ref	Size	Pairs
41455	5.5	25 x 2 / Box	41475	7.5	25 x 2 / Box
41460	6	25 x 2 / Box	41480	8	25 x 2 / Box
41465	6.5	25 x 2 / Box	41485	8.5	25 x 2 / Box
41470	7	25 x 2 / Box			

## Synthetic Surgical Gloves

### Biogel® Skinsense® Not available in New Zealand

A sterile, synthetic, polychloroprene surgical glove that has a longer chemical breakthrough time than other Biogel gloves in the presence of common surgical chemicals and cytotoxic agents\*\*.



Ref	Size	Pairs	Ref	Size	Pairs
50955	5.5	50 / Box	50975	7.5	50 / Box
50960	6	50 / Box	50980	8	50 / Box
50965	6.5	50 / Box	50985	8.5	50 / Box
50970	7	50 / Box	50990	9	40 / Box

### Biogel® Skinsense® Indicator® Underglove

Not available in New Zealand

A blue, sterile, synthetic underglove made of polychloroprene that can be used with a Biogel PI overglove to create a coloured Puncture Indication System.



Ref	Size	Pairs	Ref	Size	Pairs
40655	5.5	50 / Box	40675	7.5	50 / Box
40660	6	50 / Box	40680	8	50 / Box
40665	6.5	50 / Box	40685	8.5	50 / Box
40670	7	50 / Box	40690	9	40 / Box

Biogel®

Hands Deserve Better™

Mölnlycke®