


# Mölnlycke® Pulsed Lavage

Truly customised orthopaedic  
solutions



Battery version



AC version

Mölnlycke® Pulsed Lavage is DEHP free

Our single-use high pressure pulsed lavage system effectively cleans bone and removes debris during total joint procedures. This increases cement penetration into the cancellous bone for potentially stronger implant fixation and reduced revision rates. It is also an effective mechanical wound debridement system for cleaning wounds, removing dirt, necrotic tissue and debris.

## Available in two versions

- AC adapter: reduced environmental impact thanks to lack of alkaline batteries. Single-use handpiece and re-usable AC adapter
- Built-in battery: a flexible and mobile version with battery pack integrated into the handpiece, for quick and easy battery removal

## All included means no hidden costs

- Ergonomic and comfortable handpiece
- Short fan spray tip with double splash shield
- Long femoral canal tip
- Integrated suction and irrigation
- Cleaning rod for solving potential occlusion in the femoral tip

## Single packed or included in Mölnlycke® ProcedurePak®

- Include Pulsed Lavage in your customised procedure-specific pack for efficiency gains
- Also available as a sterile single packed

## Filtered revision Tip

- For hip revision surgery we have a special tip to handle small bone and cement parts. Available as sterile single packed

\* Pulsed Lavage's legal manufacturer is Kaiser Technology co., Ltd. Mölnlycke Health Care is importer and distributor.

  
**Mölnlycke®**

## Mölnlycke® Pulsed Lavage | AC version

### Sterile single packed

Product code	Description
9750000	Pulsed Lavage System, AC Power
9752030	AC Power Supply (AC100-240V), AUS
9751000	Coaxial Canal Brush Tip
9751001	Cone Tip
9751002	Coaxial Canal Tip
9751003	Hip Revision Tip



AC version



AC Adapter



Coaxial Canal Brush Tip



Hip Revision Tip

### Available in ProcedurePak

Product code	Description
2322130-00	Pulsed Lavage System, AC Power
2320481-00	Coaxial Canal Brush Tip



Cone Tip




Coaxial Canal Tip


## Mölnlycke® Pulsed Lavage | Built-in Battery

### Sterile single packed

Product code	Description
9750003	Pulsed Lavage System, Built-in Battery
9750004	Coaxial Canal Brush Tip
9750005	Cone Tip
9750006	Coaxial Canal Tip
9750007	Hip Revision Tip




Built-in Battery




Coaxial Canal Brush Tip

### Available in ProcedurePak


Product code	Description
2324050-00	Pulsed Lavage System, Built-in Battery
2325431-00	Coaxial Canal Brush Tip



Cone Tip



Hip Revision Tip



Coaxial Canal Tip

Find out more at [www.molnlycke.com.au](http://www.molnlycke.com.au) | [www.molnlycke.co.nz](http://www.molnlycke.co.nz)

You will be required to input the following code **1FtYrV** or your **AHPRA registration number** to access our websites.

Mölnlycke, 12 Narabang Way, Level 4, Belrose, NSW, 2085. Phone 1800 005 231. New Zealand Orders & Enquiries 0800 005 231, [www.molnlycke.co.nz](http://www.molnlycke.co.nz). The Mölnlycke and ProcedurePak trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. © 2025 Mölnlycke Health Care. All rights reserved. AU\_001816

