

Avance® Solo and Avance® Solo Adapt

Product descriptions and units of supply




Avance® Solo

– Closed incisions and open wounds (ciNPWT and dNPWT)

1 Pump, canister and border dressing

The three sets available in the Avance® Solo assortment are designed for convenience in the OR, primarily for closed incision applications. These sets contain everything required for application and commencement of therapy without duplications, so as to avoid waste. Each set contains a pump and cannister, a dressing with fixation strips and 4 x AA batteries.



2 Avance Solo pump and canister


The Avance Solo pump is always packaged in a format that includes;

1. The Pump body housing the motor and pump assembly
2. A plastic belt clip
3. A 50ml replaceable canister with tubing and female connector
4. 4 x AA Lithium batteries



3 Border dressing

The eight different sizes of dressings are supplied in retail packaging containing two individually packaged sterile dressings. Each dressing has an affixed dressing tube and male connector, and six polyurethane film, acrylic adhesive secondary fixations strips.



4 Avance Solo canister and foam


Canisters and foam wound filler are supplied in individually pouched and sterilised pieces that are combined as a group of four in a single package. Canisters have a pre-affixed tubing and connector.



Avance® Solo Adapt – Pressure injuries

1 Avance® Solo Adapt

Avance Solo Adapt couples the same Avance Solo pump with a customisable Safetac® Technology film and foam dressing. The dressing can be cut and customised to support application over difficult anatomy ie the sacral region.



2 Avance Solo Adapt Dressing

Dressings are supplied in an individual sterile pouch containing a Safetac Technology adhesive film dressing, Avance Solo Adapt Foam, and tubing with connector.



x5

3 Safetac Technology Film dressing

Safetac Technology adhesive film dressings are supplied in an individual sterile pouch.



x5



System product nomenclature

- All system products are supplied in an individual sterile pouch.
- Avance Solo and Avance Solo Adapt are supplied in a modular format to allow the user to select according to their requirements and, to minimise potential for over stocking or waste.
- Individual sterile units; products combined with a single sterile pouch or blister pack.
- Retail pack: a retail pack may contain more than one individual sterile unit.

Avance® Solo Assortment

Product code	Description	Individual Sterile Unit Contents	Pcs/Box
Pump, canister and border dressing			
8881020	Avance® Solo pump, canister and border dressing	Pump + canister 50ml with tubing + 4 x AA batteries + 10 x 20cm dressing + belt clip	1
8881030	Avance Solo pump, canister and border dressing	Pump + canister 50ml with tubing + 4 x AA batteries + 10 x 30cm dressing + belt clip	1
8881035	Avance Solo pump, canister and border dressing	Pump + canister 50ml with tubing + 4 x AA batteries + 10 x 35cm dressing	1
Avance Solo pump and canister			
8880050	Avance Solo pump and canister	Pump + canister 50ml with tubing + 4 x AA batteries + belt clip	1
Border dressing			
881020	Avance Solo border dressing – 10 x 20cm	10 x 20cm border dressing + 6 secondary fixation strips	2
881030	Avance Solo border dressing – 10 x 30cm	10 x 30cm border dressing + 6 secondary fixation strips	2
881035	Avance Solo border dressing – 10 x 35cm	10 x 35cm border dressing + 6 secondary fixation strips	2
881515	Avance Solo border dressing – 15 x 15cm	15 x 15cm border dressing + 6 secondary fixation strips	2
881520	Avance Solo border dressing – 15 x 20cm	15 x 20cm border dressing + 6 secondary fixation strips	2
881530	Avance Solo border dressing – 15 x 30cm	15 x 30cm border dressing + 6 secondary fixation strips	2
882020	Avance Solo border dressing – 20 x 20cm	20 x 20cm border dressing + 6 secondary fixation strips	2
882525	Avance Solo border dressing – 25 x 25cm	25 x 25cm border dressing + 6 secondary fixation strips	2
Other system consumables			
880050	Avance Solo canister 50ml	Canister 50ml with tubing	4
882000	Avance Solo foam 1.5 x 10 x 12cm	Foam 1.5 x 10 x 12cm	4
Avance Solo accessories			
890003	Carrier case	17 x 11 x 8cm	2

Avance Solo Adapt Assortment

Product code	Description	Individual Sterile Unit Contents	Pcs/Box
9994152	Avance Solo Adapt pump, canister, foam, film and transfer port	Pump + canister 50ml with tubing + 4 x AA batteries + Safetac® film dressing + foam + belt clip	1
994152	Avance Solo Adapt foam, film and transfer port	Safetac film dressing + foam + transfer port and tubing	5
994000	Avance Solo Adapt film	Safetac adhesive film dressing	5

Find out more at www.molnlycke.com.au | www.molnlycke.co.nz

Mölnlycke, 12 Narabang Way, Level 4, Belrose, NSW, 2085. Phone 1800 005 231. New Zealand Orders & Enquiries 0800 005 231, www.molnlycke.co.nz. The Mölnlycke, Avance and Safetac trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. © 2026 Mölnlycke Health Care AB. All rights reserved. AUWC_002288

