

Mepilex® Border Flex

Flex Technology – stays on and conforms

- 360° flexibility
- Improves comfort¹⁻³
- Minimises risk of detachment⁴
- Reduces skin stress¹⁻³



Retention layer

- Effectively manages both normal and viscous exudate⁵
- Traps exudate containing bacteria⁶
- In combination with the spreading layer and backing film, contributes to the dressing's high moisture vapour transmission rate⁵



Foam layer

- Absorbs exudate and transports it to the spreading layer^{7,8}



Highly breathable backing film

- High moisture vapour transmission rate⁷
- Exudate Progress Monitor lets you track and record exudate progress⁹
- Allows you to check colour of exudate¹⁰
- Works as a barrier against bacteria^{**11}
- Water-resistant so your patients can shower¹²



Spreading layer

- Distributes exudate over a wide surface area to maximise exudate transport to the retention layer and backing film¹³



Safetac® wound contact layer

- Minimises pain at dressing changes¹⁴ and reduces risk of maceration^{15,16}
- It's easy to handle – and you can remove the dressing without damaging the skin¹⁴
- Safetac® also protects new tissue, so wounds remain undisturbed to support natural healing^{15,16}

**particles/ microbes larger than 25 nm

Stays on. Uniquely conforms.

- Flex Technology provides excellent flexibility¹⁷ and conformability to avoid early detachment⁴
- The Exudate Progress Monitor allows the spread of exudate to be objectively assessed and recorded, without disturbing the wound⁹
- For use on a wide range of exuding wounds
- Also available in oval shape for greater clinical choice
- The 5-layer construction of Mepilex® Border Flex works to absorb, channel and trap exudate containing bacteria away from the wound bed and prevent re-entry, even under compression⁶
- Test shows that Mepilex Border Flex releases up to 22 times less bacteria than other dressings^{*18}

Easy application



The three-part release liner can help you to place the dressing accurately and securely, avoiding wrinkles, rolled edges and re-applications¹⁹.

***Available in square 15 x 15cm, 15 x 20cm and Oval 13x16 and 15x19



Mepilex® Border Flex is available in an oval shape because every wound is different, and sometimes a special shape might be just what you need for secure fixation.

Designed to conform and stay on

Mepilex® Border Flex has greater conformability compared to other bordered foam dressings tested due to unique Flex Technology¹⁷. A dressing that moves with the patient and stays on even in the shower¹⁴.

Smart exudate management

The 5-layer dressing absorbs, channels and traps exudate and allows you to track exudate progress^{6,8,9,13}. Keeping exudate away from the wound bed, and preventing re-entry, even under compression⁶.

Ordering information (sterile packed)



Square/Rectangle

Product code	Size (cm)	Pcs/box
595211	7.5 x 7.5	10
595311	10 x 10	10
595011	12.5 x 12.5	10
595411	15 x 15	10
595611	15 x 20	10



Oval

Product code	Size (cm)	Pcs/box
583500	7.8 x 10	5
583300	13 x 16	5
583400	15 x 19	5

Reduced pain and damage

Our proven Safetac® Technology means less pain and damage during dressing changes¹⁴, and reduced risk of maceration than traditional (adhesive) dressings^{15,16}.

Areas of use

Mepilex Border Flex is designed for a wide range of exuding wounds, such as pressure injuries, leg and foot ulcers, traumatic wounds (e.g. skin tears) and surgical wounds. Mepilex Border Flex can also be used on dry/necrotic wounds in combination with gels. Mepilex Border Flex reduces post-operative blistering, and may also be used as part of a prophylactic therapy to help prevent skin damage, e.g. pressure injuries.

Note: In the case of signs of clinical infection, consult a healthcare professional for adequate infection treatment.

When you scan this QR code, you will be taken directly to my.molnlycke.com.au / my.molnlycke.co.nz, where you will find:



Australia

- Product details
- Application guide and video
- Patient information guide
- Instruction leaflet



New Zealand

References: 1. Molnlycke Health Care. Data on File. PD-529747. 2018. 2. Molnlycke Health Care. Data on File. PD-534569. 2018. 3. Molnlycke Health Care. Data on file. PD-534571. 2018. 4. Molnlycke Data on file. PD-530246. 2017. 5. Molnlycke Health Care. Data on file. PD-527642. 2019. 6. Molnlycke Health Care. Data on file. PD-537072. 2019. 7. Molnlycke Health Care. Data on File. PD-524349. 2019. 8. Molnlycke Health Care. Data on file. PD-532095. 2019. 9. Molnlycke Health Care. Data on file. PD-524349. 2019. 10. Molnlycke Health Care. Data on file. 2018. 11. Molnlycke Health Care. Product Manual Mepilex® Border Flex - Viral penetration. Data on file. 12. Molnlycke Health Care. Mepilex® Border Flex - Waterproofness. Data on file. 13. Molnlycke Health Care. Data on file. PD-528874. 2018. 14. White R. A multinational survey of the assessment of pain when removing dressings. Wounds UK. 2008;4(1):14-22. 15. Woo K, Coutts PM, Price P, Harding K, Sibbald RG. A randomised crossover investigation of pain at dressing change comparing foam dressings. Adv Skin Wound Care. 2009;22(7):304-10. 16. Meaume S, et al. A study to compare a new self-adherent soft silicone dressing with a self-adherent polymer dressing in stage II pressure ulcers. Ostomy Wound Manage. 2003;49(9):44. 17. Molnlycke Health Care. Data on file. PD PD-528870. 2018. 18. Karlsson S, et al. Bacterial trapping of a newly developed all-in-one soft silicone foam dressing. Poster presented at: EWMA; 2018 May 9-11; Krakow, Poland. 19. Molnlycke Health Care. Data on file. PD-530091. 2019.

Find out more at www.molnlycke.com.au | www.molnlycke.co.nz

You will be required to input the following code **1FtYrV** or your **AHPRA registration number** to access our websites.

Molnlycke, 12 Narabang Way, Level 4, Belrose, NSW, 2085. Phone 1800 005 231. New Zealand Orders & Enquiries 0800 005 231. www.molnlycke.co.nz. The Molnlycke, Mepilex and Safetac trademarks, names and logos are registered globally to one or more of the Molnlycke Health Care Group of Companies. © 2025 Molnlycke Health Care AB. All rights reserved. AUWC_001675

