

Mölnlycke® Plume Evacuation Pencils

With an integrated telescopic function providing an immediate, on demand, solution for both deep and shallow incisions



Surgical staff are at risk every time surgical plume is created in the operating room. The solution is in your hands with the new, unique Mölnlycke® Plume Evacuation Pencils, a superior solution that ensures minimised plume exposure and clear visualisation of the operative site. Putting the safety of you and your patients first.

- A unique all in one light-weight, slim design with an integrated telescopic function providing an immediate, on demand, solution for both deep and shallow incisions.
- Mölnlycke Plume Evacuation Pencils offer a sustainable solution, DEPH and PVC free.
- Very high suction capacity (95 l/min) minimising plume exposure and providing rapid visualisation of the operative site.
- Available both in your customised surgical trays as well as single packed.

Technical data

Choice of electrodes:
Stainless steel ❶ for quick, simple cases
(only available into Procedure packs).
Coated PTFE ❷ for reduced interruption
of surgery due to eschar build up.
Insulated PTFE ❸ for reduced risk of inadvertent
tissue damage when operating in tight spaces.



360° swivel handle enables freedom of movement reducing the risk of wrist fatigue

Universal 22mm connector to fit all plume evacuation machines

Narrow diameter ensures precise control and greater visualisation, especially in restricted spaces

Convenient 4 metre cable

Ordering information (for single sterile)

Ref. No.	Product name and description
420100-00	Plume Evacuation Pencil, PTFE electrode
420101-00	Plume Evacuation Pencil, PTFE Insulated electrode

Legal manufacturer: Prima Medical Limited.
Imported and distributed by: Mölnlycke Health Care AB

Find out more at: www.molnlycke.com.au or www.molnlycke.co.nz

Mölnlycke, 12 Narabang Way, Level 4, Belrose, NSW, 2085. Phone 1800 005 231. New Zealand Orders & Enquiries 0800 005 231, www.molnlycke.co.nz.
The Mölnlycke, BARRIER and ProcedurePak trademark, name and logo are registered globally to one or more of the Mölnlycke Health Care group of companies.
©2026 Mölnlycke Health Care AB. All rights reserved. AU.002486

