

Mepilex® Lite with Safetac® Technology


One dressing, customized to fit.




reddot winner 2022

Mepilex® Lite dressings help protect fragile skin and may be used for preventing skin damage under medical devices.^{1,2} An innovative cutting guide incorporated directly on the sterile release liner of the new 5" x 5" size provides guidance on multiple cutting options to protect anatomic areas at risk of damage due to prolonged use of medical devices, such as non-invasive respiratory masks and oxygen tubing.*

New Size



Suggested uses




*The prophylactic use of dressings in reducing the risk of medical device related pressure injuries is well reported, however their use under PPE in terms of maintenance of a seal and potential impact on viral transmission has not been tested by Mölnlycke or others to our knowledge.

Mepilex® Lite


Mölnlycke®

Mepilex® Lite dressings conform well to difficult-to-dress areas³

Safetac® Technology is the original less-pain contact layer with silicone adhesion, reducing pain and risk of skin damage for patients at dressing removal.⁴



Traditional adhesive sticks only to the tops of the skin, sticking hard where it sticks. When the dressing is removed, skin may be stripped.



Safetac® technology molds completely to skin pores. At dressing removal, peel forces are distributed over a wider surface area. This minimizes the stripping.

Safetac®
TECHNOLOGY

Ordering information

Product ID	Description	Pcs/Box	Pcs/Case	HCPCS Code
284123	Mepilex Lite for use with medical devices	5	50	A6210

Further information on using Mepilex® Lite with medical devices can be accessed by **scanning the QR code on the dressing unit pack**



References: 1. Boesch RP, Myers C, Garrett T et al. Prevention of tracheostomy-related pressure ulcers in children. Pediatrics. 2012;129: e792–e797. 2. Cohen LP, Ovadia-Blechman Z, Hoffer O, Gefen A. Dressings cut to shape alleviate facial tissue loads while using an oxygen mask. Int Wound J. 2019; 1–14. 3. Mölnlycke Health Care. Data on file. 2005. 4. White, R. A multinational survey of the assessment of pain when removing dressings. Wounds UK 2008;4(1):14–22.

The Mölnlycke, Mepilex and Safetac trademarks, names and logo types are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Distributed by Mölnlycke Health Care US, LLC, Peachtree Corners, Georgia 30092. © 2022 Mölnlycke Health Care AB. All rights reserved. 1-800-843-8497. USMM_270671 MHC-2022-79877

Mölnlycke®